



BLOOD DONOR PHLEBOTOMY DEVICE EVALUATION FORM

Date: _____ Department: _____
 Occupation: _____ Product: _____ Number of times used: _____

Please **circle** the most appropriate answer for each question. Not applicable (N/A) may be used if the question does not apply to this particular product.

	Agree.....Disagree
1 The safety feature does not obstruct vision of the tip of the needle.	1 2 3 4 5 N/A
2 This product does not require more time to use than a non-safety device.	1 2 3 4 5 N/A
3 Use of this product requires you to use the safety feature.	1 2 3 4 5 N/A
4 The safety feature works well with a wide variety of hand sizes.	1 2 3 4 5 N/A
5 The device is easy to handle while wearing gloves.	1 2 3 4 5 N/A
6 There is a clear and unmistakable change (audible or visible) that occurs when the safety feature is activated.	1 2 3 4 5 N/A
7 The safety feature operates reliably.	1 2 3 4 5 N/A
8 The needle is immediately shielded or retracted upon removal from vein.	1 2 3 4 5 N/A
9 The needle shielding is engaged from behind the needle.	1 2 3 4 5 N/A
10 After being placed in the permanently locked position the safety feature cannot be undone.	1 2 3 4 5 N/A
11 This safety product is no more difficult to use than non-safety products.	1 2 3 4 5 N/A
12 The user does not need extensive training for correct use of the product.	1 2 3 4 5 N/A
13 The design of the product suggests proper use.	1 2 3 4 5 N/A
14 It is not easy to skip a crucial step in proper use of the device.	1 2 3 4 5 N/A
15 The product can be easily used in either hand.	1 2 3 4 5 N/A
16 This device provides a better alternative to traditional blood donor phlebotomy devices.	1 2 3 4 5 N/A
17 The product is compatible for use with current blood collection sets produced by a variety of manufacturers.	1 2 3 4 5 N/A

Of the above questions, which three are the most important to your safety when using this product?

Are there other questions which you feel should be asked regarding the safety features of this product?